


# A New Technology for Preparing Monoclonal Antibodies to Molecules Expressed on the Cell Surface

Paul B. Fisher

Surface-epitope masking (SEM) is a new procedure that permits the efficient and selective production of monoclonal antibodies (MAbs) reactive with both known and unknown molecules expressed on the cell surface. Applications of this technology include the development of MAbs reacting with the multidrug resistance gene, human interferon receptors, novel human prostate carcinoma tumor-associated antigens (TAAs), and human breast carcinoma TAAs. In theory, the SEM approach can be used to develop MAbs and to identify genes associated with important cellular processes including growth, differentiation, immunological recognition, tumorigenesis, metastasis, cell senescence, atypical MDR, and autoimmune diseases.

**Paul B. Fisher, PhD**, is a professor of clinical pathology and a Chernow research scientist in the Departments of Pathology and Urology at the College of Physicians and Surgeons, Columbia University, 630 West 168th Street, New York, NY 10032. He is also the director of neuro-oncology research at Columbia-Presbyterian Medical Center, New York, NY, tel. (212) 305-3966 or 305-7870, fax (212) 305-8177. **Sidney Pestka, MD**, column editor, is professor and chair of the Department of Molecular Genetics and Microbiology, UMDNJ-Robert Wood Johnson Medical School, 675 Hoes Lane, Piscataway, NJ 08854-5635.



The production of monoclonal antibodies (MAbs) targeted toward defined cell surface-expressed molecules is often an inefficient, time-consuming, and expensive procedure (1-3). Recently, a method based on immunological subtraction was devised to circumvent these problems. This approach is called surface-epitope masking (SEM) (4).

The SEM process involves the generation of polyclonal antibodies reacting with a complete and random set of surface antigens expressed on one cell type (referred to as a *driver*). These polyclonal antibodies coat (mask) common surface antigens on a related cell type expressing a surface molecule to be targeted for the production of MAbs (referred to as a *tester*). The SEM procedure allows efficient generation of hybridomas that produce MAbs reacting with cell surface-expressed molecules of both known and unknown functions (4).

A major problem that often limits the effectiveness of chemotherapy is the acquisition by tumor cells of cross-resistance to structurally diverse chemotherapeutic agents, a process called *multidrug resistance* (MDR) (5,6). The genetic lesion that induces one form of MDR — namely, typical MDR — involves increased expression of a gene (*mdr-1*) that encodes a 170,000 molecular weight cell surface-expressed drug-efflux transporter protein (P-glycoprotein) (5).

Studies have determined whether the SEM strategy could develop MAbs that react specifically with external surface-localized epitopes of the human MDR P-glycoprotein (4). In this process, the driver cell is a cloned rat embryo fibroblast cell line, CREF-Trans 6 (6,7), and the tester cell line is a genetically altered MDR-CREF-Trans 6 cell clone. MDR-CREF-Trans 6 cells are incubated with high-titer

polyclonal antibodies produced against the untransfected CREF-Trans 6 parental cells. Surface-epitope-masked MDR-CREF-Trans 6 cells are injected into BALB/c mice, and immune spleen cells taken from these mice are fused with myeloma cells. This process results in the development of hybridomas that produce MAbs that react by fluorescence activated cell sorter (FACS) analysis with tester MDR-CREF-Trans 6 cells, but not with driver non-MDR-CREF-Trans 6 parental cells (4). The ability of these genetically altered MDR-CREF-Trans 6 cells to react with the SEM-derived MDR MAbs by FACS indicates that these antibodies can interact with specific surface-expressed native epitopes on the MDR P-glycoprotein.

To determine the utility and reactivity of the SEM-derived MDR MAbs toward additional cell types expressing the human typical MDR gene, the human breast carcinoma cell line MCF-7 was genetically modified with the human *mdr-1* gene (4). Whereas the non-MDR parental MCF-7 cells do not react with the SEM-derived MDR MAbs, all of the independently derived MDR-MCF-7 clones are reactive by FACS analysis with the SEM-derived MAbs. These findings indicate that the SEM approach can target MAb production toward epitopes on a defined human surface-expressed protein.

A further application of the SEM approach for developing MAbs that react with known surface-expressed proteins involves human interferon receptors. In a collaborative project with Dr. Sidney Pestka's laboratory (Robert Wood Johnson Medical School, University of Medicine and Dentistry of New Jersey), CREF-Trans 6 cells were genetically engineered to express either human IFN- $\gamma$  or IFN- $\alpha$  receptors. CREF-Trans 6 cells expressing the IFN- $\gamma$  receptors were used as tester cells for application of the SEM procedure. Hybridomas were produced that secrete MAbs reacting with human IFN- $\gamma$  surface receptors in genetically engineered CREF-Trans 6 cells, human cells, and somatic cell hybrids containing the appropriate human genes and encoding the human IFN- $\gamma$  receptor. A similar SEM approach is currently being used to develop MAbs that react with human IFN- $\alpha$  receptors and surface-expressed human interferon accessory proteins. These findings provide further scientific support for SEM and for the use of this approach to tailor the development of MAbs that react with specific surface-expressed molecules.

The SEM approach also has wide applicability to the development of MAbs that react with surface-expressed molecules with unknown functions. For example, with human prostate carcinoma DNA transfected nude mouse tumor-derived CREF-Trans 6 cells, the

SEM procedure produced MABs that react with surface molecules expressed on transfected (genetically altered) but not on untransfected CREF-Trans 6 cells (4). These MABs, known as Pro 1.1 to Pro 1.5, also react with the surface of human prostate carcinoma cell lines, including LNCaP (initial source of transforming DNA), DU-145, and PC-3 (4). In addition, Pro 1.5 reacts with surface epitopes on human breast carcinoma cells, including MCF-7 and T47D, but not with human skin fibroblasts, specific human colon carcinomas, or melanomas (4). In recent studies, Pro 1.5 has displayed reactivity with fresh-frozen sections of patient-derived prostate carcinomas, but not with normal prostate or benign prostatic hyperplasia tissue.

A similar approach was used to develop MABs that react with the surface of human breast carcinomas. High molecular weight DNA from the human breast carcinoma cell line T47D was cotransfected with an antibiotic resistance gene (pSV2neo) into CREF-Trans 6, and G418-resistant colonies were injected into athymic mice. Tumors that developed were removed, established in culture, and subjected to the SEM protocol. The

SEM-derived MABs, Br-car 4.2.1 and Br-car 5.2.4, displayed reactivity with the surface of appropriately transfected CREF-Trans 6 cells, as well as T47D and MCF-7 human breast carcinoma cells. Br-car 4.2.1 and Br-car 5.2.4 selectively detected breast carcinomas in both fresh-frozen and formalin-fixed patient-derived specimens. These studies indicate the power of the SEM approach in developing immunological reagents with direct diagnostic potential.


The SEM-derived MABs can also be used to clone genes that encode unidentified surface-localized molecules expressed by genetically modified CREF-Trans 6 cells. For example, screening of an LNCaP cDNA expression library with Pro 1.5 resulted in cloning of prostate carcinoma tumor-antigen gene-1 (PCTA-1) (unpublished data). PCTA-1 is a novel gene, not reported in current DNA databases, of 3.8 Kb that encodes a protein with a predicted  $M_r$  of ~36 Kd. A protein of similar size is detected after in vitro translation of the PCTA-1 cDNA in a rabbit reticulocyte translation system and by immunoprecipitation from  $^{35}$ S-methionine-labeled lysates from LNCaP and DU-145 cells. With specific sequences present in PCTA-1, primers are

now available that detect PCTA-1 expression by RT-PCR in transfected cells, prostate carcinoma cell lines, and patient-derived prostate carcinomas. These results indicate that SEM can develop MABs of clinical interest and can identify genes that encode these proteins.

The SEM protocol may be useful for enhancing the sensitivity of the "phage display combinatorial library" approach to identifying genes involved in transformation and genes that encode surface-expressed human tumor-associated antigens (8,9). For example, human tumor DNA-transfected CREF-Trans 6 polyclonal antibody and used as immunogens to stimulate an immune response in mice. Isolated spleen cells can then directly generate combinatorial phage cDNA libraries. In this way, genes that encode surface-expressed molecules involved in transformation and cell surface-expressed, tumor-associated antigens can be directly identified and cloned.

#### POTENTIAL APPLICATIONS


Conceptually, SEM represents a broad-based "immunological subtraction" technology. It will permit the directed development of MABs and the identification of genes reacting with




**Reliability.  
Compliance.  
Calibration Manager.**

*Calibration Manager™* is the leading calibration management software designed specifically for companies regulated by the FDA. Over 800 installations world-wide rely on *Calibration Manager* to control and maintain their critical calibration records. If you want to go through your next inspection with confidence, use *Calibration Manager*.

Blue Mountain Software, Inc.  
208 W. Hamilton Avenue  
State College, PA 16801  
Fax: (814) 234-7077 **Circle 3**



**For a  
FREE  
Demo Disk  
Call:  
800-  
982-2388**



*Over 20 Years of  
Custom Pharmaceutical  
Manufacturing Service*

**BELL-MORE LABS, INC.**  
Freeze Dried Parenterals Our Speciality  
4030 Gill Avenue  
Hampstead, MD 21074 U.S.A.  
410-239-7554

and encoding, respectively, diverse cell surface-expressed molecules present in a tester target cell but missing in a driver target cell. Although our current use of SEM has involved the CREF-Trans 6 system, the applications of this approach are much broader. Numerous additional applications of the SEM approach are possible. Following are some examples of these applications.

**Developing MABs specific for early-stage and metastatic cancers.** Polyclonal antibodies

are produced against an early-stage or metastatic tumor of the same histology, followed by the SEM approach to develop, respectively, early-stage- or metastatic-specific MABs. These SEM-derived MABs would be useful for diagnostics or, if they displayed appropriate specificity, for therapeutic applications.

**Identifying genes encoding surface-expressed molecules found preferentially in early-stage and metastatic cancers.** The SEM approach can be coupled with the com-

binatorial phage cDNA expression library approach to identify cDNAs that encode surface-expressed early-stage cancer or metastatic cancer antigens. The cDNAs can be sequenced and this information can be used to develop primers for RT-PCR, thereby resulting in potentially important cancer diagnostic tools.

**Developing MABs reacting with patient-specific, tumor-associated, cell surface-expressed molecules.** To achieve this goal, polyclonal antibodies produced against adjacent normal tissue (e.g., normal breast tissue) mask antigens on genetically matched tumor tissue (e.g., breast carcinoma). SEM-derived MABs reacting with specific tumor-associated antigens on a patient's carcinomas (which may also represent generic antigens expressed on additional patient breast carcinomas) may prove useful for identifying metastatic spread and for therapeutic applications (immunoconjugates with radio-nuclides or toxic drugs).

**Identifying genes encoding patient-specific, tumor-associated, surface-expressed molecules.** The SEM approach described in the previous paragraph is combined with the combinatorial phage cDNA expression library approach to identify cDNAs that encode surface-expressed, patient-specific, tumor-associated antigens.

**Targeted development of MABs specific for the outer domain of membrane-localized growth factor receptors and cell-membrane transporter proteins.** In this application, polyclonal antibodies produced against a nonexpressing parental cell line (driver) mask the same target cell genetically altered to express the surface membrane protein (tester), e.g., bombesin receptor, interferon receptors, and other proteins. SEM-derived MABs would prove useful for both diagnostic and functional studies.

**Development of primate monoclonal antibodies.** In this approach, genetically altered CREF-Trans 6 cells are coated with high-titer CREF-Trans 6 polyclonal antibodies and injected directly into primates. Primate MABs developed from immune-stimulated primate spleen cells should prove less immunogenic in humans than murine monoclonal antibodies (3).

#### CONCLUSION


With appropriate modifications, the SEM approach should also prove useful for producing MABs and/or isolating relevant genes involved in tumor cell recognition by both nonspecific and specific immunologic effector cells, reacting with and encoding novel native and chimeric tumor-associated growth factor receptors, reacting with and encoding cell surface-localized tumor suppressor gene

## OUR "QC" SERIES ULTRA-SANITARY butterfly valve...



**offers more, but doesn't cost more!**

- ☐ more styles
- ☐ more sizes
- ☐ more designed-in features
- ☐ more accessories
- ☐ more mounting methods
- ☐ more quick-connect types



EASIER TO DISASSEMBLE AND CLEAN.

Plus 20 years experience working with high potency, sterile and aseptic products.

**Yes we have competitors, but not competition!**

WRITE FOR LITERATURE


# CORA

INTERNATIONAL • INC

476 ROUTE 9, ENGLISHTOWN, NEW JERSEY 07726  
(908) 536-6600 • FAX (908) 972-2630

599-504

# STERILIZERS, DRYERS, AND CONTINUOUS PROCESS



Gruenberg's 60-year history of building quality dry-heat sterilizers and granulation dryers provides a firm foundation for new directions.

Sterilizers themselves can be customized to include automatic doors, pressure balancing, and PLC control systems. Both sterilizers and dryers can be integrated with hardwall insulation or laminar flow systems.

If you need combined expertise in temperature, cleanliness, airflow, and humidity control, call or write Gruenberg today.

◀ Sterilizers and dryers are available in bench-top, cabinet, truck-in, pass-through, and continuous arrangements.



2121 Reach Road, Williamsport, PA 17701  
717-326-1755 Fax: 717-326-7304

Circle 81

## Let Shaklee Manufacturing Excellence Work For You

Extend your operational capabilities with Shaklee Contract Manufacturing. We'll customize our services to deliver solid dosage pharmaceuticals to fit your most exacting requirements—and your budget.

Our experienced staff and qualified facilities are in place and ready to assist every step of the way. We'll handle development, testing, technology transfer, validation, manufacturing and packaging.

Shaklee is known for quality, reliability, and personal service. We understand your need to reduce time-to-market. We offer high volume capability and an on-time delivery record that is second to none.

For more information on how we can work for you, please contact Tim Reason at (510) 887-5179.



Part of the  
Yamanouchi  
group of companies.

*Our pilot plant is maintained to cGMP standards for all phases of clinical supply manufacturing.*

Circle 60

products, mediating atypical multidrug resistance, and identifying mediators of autoimmune diseases.

Overall, SEM is a powerful new method of preparing monoclonal antibodies to tumor-associated antigens and other surface proteins of cells.

### ACKNOWLEDGMENT

The author gratefully acknowledges the contributions of Drs. Neil I. Goldstein, Jorge Leon, Jian Lin, Sidney Pestka, Ruqian Shen, and Zao-zhong Su to the studies briefly reviewed in this article. Research support was provided in part by Public Health Service grants CA35675 and CA43208 from the National Cancer Institute and by the Chernow Endowment.

### REFERENCES

1. T.A. Waldman, "Monoclonal Antibodies in Diagnosis and Therapy," *Science* 252, 1657-1662 (1991).
2. L.A. Pirofski et al., "Current State of Hybridoma Technology," *Am. Soc. Microbiol. News* 58, 613-617 (1992).
3. J.A. Leon et al., "New Approaches for the Development and Application of Monoclonal Antibodies for the Diagnosis and Therapy of Human Cancer," *Pharmacol. Therapeut.* 61, 237-278 (1994).
4. R. Shen et al., "Surface-Epitope Masking: A Strategy for the Development of Monoclonal Antibodies Specific for Molecules Expressed on the Cell Surface," *J. Natl. Cancer Inst.* 86, 91-98 (1994).
5. M.M. Gottesman and I. Pastan, "Biochemistry of Multidrug Resistance Mediated by the Multidrug Transporter," *Annu. Rev. Biochem.* 62, 385-427 (1993).
6. Z.-z. Su et al., "Transfer of a Dominant-Acting Tumor-Inducing Oncogene from Human Prostatic Carcinoma Cells to Cloned Rat Embryo Fibroblast Cells by DNA-Transfection," *Anticancer Res.* 12, 297-304 (1992).
7. R. Shen et al., "Identification of the Human Prostate Carcinoma Oncogene PTI-1 by Rapid Expression Cloning and Differential RNA Display," *Proc. Natl. Acad. Sci. USA* 92, 6778-6782 (1995).
8. J.D. Marks et al., "Bypassing Immunization: Human Antibodies from V-gene Libraries Displayed on Phage," *J. Mol. Biol.* 222, 581-597 (1991).
9. W.D. Huse et al., "Generation of a Large Combinatorial Library of the Immunoglobulin Repertoire in Phage Lambda," *Science* 249, 1275-1281 (1989). ■